

interscience


Scan® 50 is a new generation colony counter, offering contrasted colony lighting while protecting the user's eyes. It improves the ergonomics of colony counting and prevents MusculoSkeletal disorders (MSD). Traceability is ensured by USB port.




Contrasting colony lighting: Dark Field 2


Dark Field 2 is an indirect lighting system with adjustable intensity. Light reflectors diffuse the light at all points without any reflections or shadows. It makes all colonies stand out on the surface, in the mass and around the edge of the dish.

Flashes of light are eliminated and the **operator's eyes are protected.**









Ergonomics redesigned

Audio Comfort

 9 different sounds and headphones jack to work without disturbing colleagues

Light signal

Light signal for each counted colony

Optimal working posture

- Inclined reading surface with magnifying lens⁽²⁾, hand rest and HandPad⁽¹⁾
- Prevents Muscular-Skeletal Disorders



Ambidextrous

- Duplicated count buttons on both sides
- HandPad, 180° arm-rest that turns on the machine

Technical specifications & Accessories


Delivered with: 1 power cord, 1 USB cable, 1 black marker pen, 1 white diffusing background,1 user manual, 1 headphones (only on Scan® 50 Pro). Magnifying lens in option.



Magnifying lens

Magnifying lens + 1 flexible arm for Scan® 50 and Scan® 50 pro

Ref. 435 001



Wolffhuegel counting grids

Wolffhuegel[™] transparent and white counting grids for Scan® 50 and Scan® 50 pro

Ref. 435 025



Spiral® counting grids (90 mm)

easySpiral® counting grids: 1 transparent and 1 white diffusing background for Scan® 50 and Scan® 50 pro - Ø 90 mm

Ref. 435 036



Circle counting grids (90 mm)

easySpiral® Dilute circle counting grids: Ø 90 mm - 1 transparent and 1 white diffusing background for Scan® 50 and Scan® 50 pro

Ref. 435 038



Adaptor for Petrifilm[™]/MC-Media Pads[™]/ Compact Dry[™]

Adaptor for Petrifilm[™], MC-media Pads[™], Compact Dry[™] for Scan[®] 50 and 50 pro

Ref. 435 051



Headphones

Headphones for Scan® 50 and Scan® 50 Pro. Connection with a jack socket.

Ref. 435 052



Optical cleaner + wipes

Optical cleaner + wipes

Ref. 435 702

Bar-code reader Bar-code reader (1D)

Ref. 522 000

Your local distributor



Scan® 50 Scan® 50 Pro Reference 435 050 435 055 Stainless steel LIMS/SIL connection USB connection Counting on pour, surface, Spiral® and circle plated Petri dishes Counting on chromogenic Petri dishes Counting on Petrifilm™, Compact Dry™, MC-Media Pads™, easyPlate™, filtration membranes LED lighting without reflections HandPad: ergonomic arm-rest 1 Headphones included Manual (0-1999 CFU) Type of colony counter Manual control to add or substract colonies Counting light signal Counting sound signals 9 different sounds / 5 volume levels Touch surface with adjustable Counting tone sensitivity Dark Field 2: Indirect lighting Lighting technology with white LED Lighting system LED lighting with 5 intensity levels SPECIFICATIONS Petri dish dimensions From Ø 55 to Ø 90 mm Via USB, LIMS, Data export CSV or text file Voltage - Frequencies 100-240 V~ 50-60 Hz 3-vear Warranty (after registration of the warranty card) In compliance with ISO 7218 and FDA BAM 27.3 x 29.5 x 27.3 x 35.2 x Dimensions (w x d x h) 15.1 cm 15.1 cm Net weight 4.05 kg 4.45 kg Box dimensions (w x d x h) 35 x 45 x 40 cm 35 x 45 x 40 cm



Gross weight







Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.

Contact us for IQOQ and warranty extension

interscience

PARIS

Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770 - Email: sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

SHANGHA

Phone: +86 (0)21-64739390 - Email: sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232 - Email: sales.asia@interscience.com

Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

