



170 Shields Court Unit 2 Markham, ON L3R 9T5 <u>TEL</u>: (905) 475-5880 ext. 226 <u>FAX</u>: (905) 475-1231





Closed Door System Fully Automatic Autoclaves



Autoclaves will meet all your sterilization needs

- Chamber door remains closed during entire cycle
- Hepa filtered air is circulated in chamber during drying
- All program parameters can be customized and stored
- New brighter and easier to read panel
- Two year parts and labor warrant
- * Above photo shown with optional printer



With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer fully automatic autoclaves fill, sterilize, exhaust and dry at the touch of a button.

Item No.	Chamber Dimensions Ø x L In(mm)	Volume	Overall Dimensions DxWxH In(mm)	Standard Cassettes	Tray Dimensions DxWxH In(mm)	No. of Trays	Standard Unwrapped Cycle time	Voltage(V) Freq.(Hz)	Power(W) Current(A)	Shipping Weight Ibs(kg)
EZ9™	9″ x 18″ (230 x 470)	5gal/19L	21.5″x20″x14.4″ (545x510x365)	2 full/2 half	16.3″x6.7″x0.8″ (415x170x20)	3	Cold: 27 mins Hot: 13 mins	120V 50/60	1400W 12A	95 (43)
EZ10™	10″ x 19″ (254 x 475)	6gal/23L	21.5″x20″x14.4″ (545x510x365)	3 full/3 half	16.3″x6.7″x0.8″ (415x170x20)	4	Cold: 30 mins Hot: 14 mins	120V 50/60	1400W 12A	95 (43)
*EZ10k™ Kwiklave™	10″ x 19″ (254 x 475)	6gal/23L	21.5″x20″x14.4″ (545x510x365)	3 full/3 half	16.3″x6.7″x0.8″ (415x170x20)	4	Cold: 21 mins Hot: 11 mins	230V 50/60	2200W 10A	95 (43)

E≈9" and E≈10" Series:

*Model EZ10k™ requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5KVA) may be required.

Features:

- · Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- High efficiency air pump circulates air in chamber during drying cycle
- 0.2µm Hepa air filter (replaceable)
- New brighter and easier to read panel
- · Program parameters can be customized and stored


Control Panel:



Optional Internal Printer



Documents: date, time, temperature, pressure.





Steam Distiller

Chamber Brite Autoclave Cleaner

Tuttnauer USA Co., Ltd., 25 Power Drive Hauppauge, NY, 11788 Tel: (800) 624 5836, (631) 737 4850, Fax: (631) 737 0720 E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

Cleaners and Tablets

