



Cat. No.	Ordering Specifications
5002	AUTOMATIC HIGH SENSITIVITY OSMETTE A™ Laboratory Model Osmometer for both 0.2 & 2.0 ml samples. Complete with instruction manual, 24 2.0 ml sample tubes, 12 0.2 ml sample tubes, spare stirring wires, 1 bottle each of 100, 500 & 1500 mOsm/kg standards and bath liquid, test tube rack. Shipped assembled and ready to use. (115V*, 50/60 Hz)

ACCESSORIES

2012	High Isolation Probe for 0.2 & 2.0 ml samples
2021	Glass Sample Tubes, 2.0 ml, 12/pk
2022	Glass Sample Tubes, 0.2 ml, 12/pk
2028	<i>CLEARVIEW Disposable Sample Tubes, 0.2ml, 500/pk</i>
2032	Stirring Wires for 2.0 & 0.2 ml tubes, 4/pk
2041	Rack for 0.2 & 2.0 ml sample tubes
2093	Instruction Manual for OSMETTE A Model 5002
2100	Bath Liquid, ready to use, 250 ml
2101	Standard Solution, 100 milliosmoles, 125 ml bottle
2103	Standard Solution, 300 milliosmoles, 125 ml bottle
2105	Standard Solution, 500 milliosmoles, 125 ml bottle
2109	Standard Solution, 900 milliosmoles, 125 ml bottle
2110	Standard Solution, 1000 milliosmoles, 125 ml bottle
2111	Standard Solution, 1100 milliosmoles, 125 ml bottle
2115	Standard Solution, 1500 milliosmoles, 125 ml bottle
2201	CON-TROL 100, 100 mOsm Reference Solution, 5 ml ampules 12/pk
2202	CON-TROL 290, 290 mOsm Reference Solution, 5 ml ampules including 10 of 290 mOsm plus 1 each of 100 & 500 mOsm
2205	CON-TROL 500, 500 mOsm Reference Solution, 5 ml ampules 12/pk
2215	CON-TROL 1500, 1500 mOsm Reference Solution, 5 ml ampules, 12/pk
2263	CON-TROL Cal Pack, 4 each: 100, 500 & 1500 mOsm - 5 ml ampules